

# Evaluating CO<sub>2</sub> Tracer Gas Ingress for Integrity Assessment of Flexible Film Packaging



## Introduction

Ensuring Container Closure Integrity (CCI) in parentals is critical for maintaining sterility and product stability in pharmaceuticals.<sup>1</sup> Flexible film packaging systems (e.g., bags and pouches) present specific challenges for Container Closure Integrity Testing (CCIT) due to material permeability, seal variability, and susceptibility to micro-defects.<sup>2 3</sup> Several CCIT technologies are currently applied to such systems, including probabilistic methods such as dye ingress and bubble emission testing, as well as deterministic techniques like High Voltage Leak Detection (HVLD) and vacuum decay. While these methods are well established, their applicability to flexible film formats may be influenced by package compliance, material permeability, product conductivity, and test sensitivity limitations for very small defects.

In this study we present how laser-based Headspace Analysis (HSA), with carbon dioxide (CO<sub>2</sub>) pretreatment can be applied to packages of flexible film. This concept with detection of ingress of CO<sub>2</sub> as a trace gas has been demonstrated for other packaging systems, like vials and is today presently used in the industry but has been limited to rigid containers. Tunable Diode Laser Absorption Spectroscopy (TDLAS) enables highly sensitive and specific detection of gases such as CO<sub>2</sub>. When CO<sub>2</sub> is used as a tracer gas, CCIT can be evaluated by measuring controlled ingress into the package headspace through defects in the film or seals.

This study evaluates CO<sub>2</sub> tracer gas ingress combined with laser-based HSA as a deterministic approach for CCI assessment of flexible film packaging systems.

### Honourable Mention

This research was conducted in collaboration with Oxford Lasers, whose expertise supported the evaluation of laser based CO<sub>2</sub> tracer gas analysis as a deterministic solution for flexible packaging integrity.




GPX1500 Film Pharma - Instrument used in study

## Method

Flexible film bags were prepared with controlled artificial defects created by precision laser drilling and certified by Oxford Lasers. Three nominal defect sizes were introduced: 10 µm, 20 µm, and 50 µm.

- All bag samples were placed in a pressure vessel and exposed to 100% CO<sub>2</sub> at 1 bar overpressure for 45 minutes to allow tracer gas ingress through the defects.
- A subset of samples was additionally subjected to 24 hours of overpressurization to evaluate permeation behavior over time.
- Internal headspace CO<sub>2</sub> concentration was measured using the Gasporox GPX1500 Film Pharma instrument using laser-based HSA technology, with measurements performed at baseline (pre-exposure) and after 1, 5, 20, and 45 minutes of CO<sub>2</sub> exposure.
- CO<sub>2</sub> concentration was reported as percent headspace concentration and evaluated as a function of defect size and exposure time.

CO<sub>2</sub> Concentration (%) as a function of CO<sub>2</sub> exposure time and defect size



## Results

The CO<sub>2</sub> ingress study demonstrates a clear correlation between defect size and gas accumulation over time under 1 bar overpressure. Bags containing a 50 µm laser-drilled defect showed rapid and substantial CO<sub>2</sub> ingress, reaching 4.8% concentration at 45 minutes. In contrast, bags with a 20 µm defect resulted in significantly lower accumulation, reaching approximately 0.22% at 45 minutes. For bags with a 10 µm defect, CO<sub>2</sub> levels increased more gradually, approaching 0.19% at 45 minutes.

Extended pressure pretreatment time demonstrated measurable contributions from film permeation, reaching 5.48% at 24 hours; however, defect-driven ingress remained clearly distinguishable from background permeation. Notably, the permeation-driven CO<sub>2</sub> ingress exhibited a non-linear time dependence, consistent with diffusion-controlled transport rather than linear accumulation over time.

Defect size	0 min	1 min	5 min	20 min	45 min	24 h
10 µm	0.00	0.00	0.01 - 0.02	0.06 - 0.07	0.18 - 0.2	
20 µm	0.00	0.01- 0.03	0.01 - 0.04	0.07 - 0.09	0.22 - 0.23	
50 µm	0.00	0.09 - 0.10	0.31 - 0.35	1.25 - 1.87	4.31 - 5.29	
No defect (permeation)	0.00	0.00	0.00	0.04 - 0.05	0.11 - 0.12	5.34 - 5.62

## Conclusion

Measurement of CO<sub>2</sub> tracer gas ingress using laser-based HSA provides a sensitive, quantitative, and deterministic method for CCIT of flexible-film packaging systems. Under controlled overpressure conditions, CO<sub>2</sub> accumulation in the package headspace correlated directly with laser-drilled defect size, demonstrating clear distinction between 10 µm, 20 µm, and 50 µm defects. The method enabled objective detection of micro-defects while distinguishing ingress from material permeation, which is critical for polymeric and semi-permeable systems. These results demonstrate that laser-based CO<sub>2</sub> tracer gas analysis is a robust, deterministic, and scalable solution for CCI testing of flexible packaging systems, offering improved sensitivity and objective measurement.

## References

- [1] United States Pharmacopeia (2017). General Chapter, (1207) Package Integrity Evaluation—Sterile Products. USP–NF. Rockville, MD: United States Pharmacopeia. [https://doi.org/10.31003/USPNF\\_M99926\\_01\\_0](https://doi.org/10.31003/USPNF_M99926_01_0)
- [2] I. Ilhan, D. Turan, I. Gibson, and R. ten Klooster, “Understanding the factors affecting the seal integrity in heat sealed flexible food packages: A review,” *Packaging Technology and Science*, vol. 34, no. 6, pp. 321–337, 2021, doi:10.1002/pts.2564.
- [3] M. Giacinti Baschetti and M. Minelli, “Test methods for the characterization of gas and vapor permeability in polymers for food packaging application: A review,” *Polymer Testing*, vol. 89, p. 106606, 2020, doi:10.1016/j.polymertesting.2020.106606

## Contact us

Gasporox AB (Publ.)  
Tellusgatan 13  
224 57, Lund  
Sweden

**Office:** +46 46 540 50 40

**Email:** [info@gasporox.se](mailto:info@gasporox.se)

**[www.gasporox.se](http://www.gasporox.se)**